

# Quantitative BOLD modeling of brain oxygenation during vasodilation

ISMRM-Endorsed Workshop – Quantitative Imaging

Linh Le

UC Davis

# Oxygen Extraction Fraction


OEF is a measurement of oxygen consumed by metabolism

- Function of perfusion and oxygen metabolism
- Can be combined with perfusion to measure oxygen metabolism
- Arterial-venous difference in blood oxygen saturation

OEF is dimensionless

- Healthy resting brain range: 0.3-0.4

$^{15}\text{O}$  PET OEF measurements



Bremmer JP, et al. Mol. Imaging Biol. 2011; 13:759-768

# Oxygen Extraction Fraction

## OEF consumption in brain


- Altered during disease/activity

## Current benchmark

- Triple oxygen PET
- Highly specialized, invasive, expensive and difficult to perform

## Alternative quantitative technique


- Clinical applicable



# Quantitative BOLD (qBOLD)

## Quantitative BOLD

- Models MR signal decay in microvessels
- qBOLD signal is influenced by concentration of deoxyhemoglobin, which is relatively sensitive to OEF
- Quantified through the reversible transverse relaxation rate  $R2'$



Cherukara et al., Neuroimage 2019; 202: 116106


# Methods

## ASE data: Asymmetric Spin Echo

- Acquires only a single TE
- Refocusing pulse is moved

## Participants

- 3 healthy controls
- Received ASE before and after acetazolamide (dose of 15mg/kg) (i.e. vasodilator)



# How do we quantify OEF?

## 1. Acquire R2'-weighted data

- Asymmetric Spin Echo
- Achieve gradient echo sampling of spin echo

## 2. Estimate R2' from long tau data


- Log-linear fit to  $\tau > 15\text{ms}$  data
- $R2' \propto \text{OEF} \times \text{DBV}$

## 3. Estimate DBV from spin echo

- Difference between intercept and measured SE signal

## 4. Estimate OEF from R2' and DBV


- Known constants of proportionality used to quantify OEF




Stone and Blockley, NeuroImage (2017); 147:79-88



# ASE Images during vasodilation



# Results - OEF decreased during vasodilation






# Results – Whole brain analysis


## Subject #1

%OEF change = **9.4%**



## Subject #2

%OEF change = **21.2%**




# Conclusion

- qBOLD is efficient in computational time, non-invasive
- There was a downward trend of OEF with increasing acetazolamide response

## Future steps:

- Linear model is noisy
  - Consider Bayesian framework

# Acknowledgement



**Thanks to**

Prof. Nick Blockley (University of Nottingham)

**UC Davis Fan's Lab**

Corinne Alison Donnay (UC Davis)

Emily Holy (UC Davis)

Greg Wheeler (UC Davis)